

# Molecular analysis of cultured CTC (circulating tumor cells) from lung cancer patients

Eunjoo Hwang<sup>1</sup>, Myoung Shin Kim<sup>1</sup>, Ji-Hyun Uh<sup>1</sup>, Dong Hyoung Lee<sup>1</sup>,  
Se-Hoon Lee<sup>2</sup>, Joo Kyung Park<sup>2</sup>, Byung Hee Jeon<sup>1</sup>, Jinseon Lee<sup>1</sup>

<sup>1</sup> Cytogen Inc., Seoul, Korea,

<sup>2</sup> Department of Internal Medicine, Seoul National University Hospital,  
Seoul National University College of Medicine, Seoul, Korea

## Abstract

Circulating tumor cells (CTCs) are cells that have shed into the bloodstream from primary tumor and thus rare. CTCs could provide better understanding of tumor metastasis and noninvasive monitoring of the disease progression. However the isolation and characterization has been a major technological challenge due to their rareness. Here, we suggest the CTC culture as an effective method to obtain ample number of CTCs, which can be used for molecular analysis of original tumor characteristics and further applications. We isolated and successfully cultured the CTCs from four lung cancer patients. We analyzed the cultured CTCs to detect preexisting ALK (anaplastic lymphoma kinase) fusion using real-time PCR method, and confirmed that cultured CTCs have maintained the molecular characteristics similar to those found in primary tumors. These results suggest that the isolation and culture of CTCs can be a substitute method for tumor tissue biopsy, and may provide clinical applications, including serial blood samplings for the personalized cancer therapy based on their genomic information.

## Introduction

Circulating tumor cells (CTCs) are present in the blood of cancer patients at low concentrations. It has been proposed that CTCs may be a prospective prognostic marker for cancer progression and survival in several types of cancer [1, 2] and a potential source of the metastatic tumor cells [3, 4]. Viable CTCs isolated from cancer patients can be a useful tool for identifying molecular targets and developing new cancer treatments. However, the isolation and characterization of CTCs are technically challenging due to their rareness and heterogeneity [5]. Lung adenocarcinoma is the most common subtype of lung cancer today. Recently, the treatment paradigm for advanced non-small cell lung cancer (NSCLC) has been transformed from conventional chemotherapy to targeted therapy based on molecular aberrations in primary tumor [6]. Now, it has been regarded as standard procedure to test lung carcinoma for the presence of *EGFR* mutation and *ALK* rearrangement upon diagnosis, in order to select patients for initial *EGFR* tyrosine kinase inhibitor (TKI) and *ALK* inhibitor therapy. However, the detection of such molecular abnormalities is complicated due to the difficulty in obtaining tumor material from repeated tissue biopsies [7]. Here, we were able to obtain sufficient amounts of CTCs through CTC cultures, and performed molecular analyses using cultured CTCs to confirm preexisting *ALK* rearrangement.

## Materials and Methods

### Blood collection


Blood samples (5-10 ml) from advanced NSCLC cancer patients

### Primary culture of CTCs

Cytogen protocol [9]

Culturing in growth medium (RPMI-1640 medium supplemented with 10% FBS, 2% antibiotic-antimycotic) at 37°C, 5% CO<sub>2</sub> for 16-18 days

### Cytogen Protocol



### Immunofluorescence analysis

EpCAM (Cell Signaling), CD45 (Santa Cruz), DAPI staining

### Immunocytochemistry

EpCAM (Cell Signaling)


### FISH (Fluorescent in situ hybridization) for ALK rearrangement

### Quantitative real-time PCR for EML4-ALK fusion detection


AmoyDx EML4-ALK Fusion Gene Diagnostic Kit (Amoy Diagnostics Company Ltd.)

## Results


Representative images of CTC culture at day 0, 4, 7, and 16 (X400).




Real-time PCR analysis of cultured CTCs for EML4-ALK rearrangement.




Immunofluorescent staining for EpCAM (green), CD45 (red), and nuclei (blue) (X400).



H&E staining (a) and Immunocytochemical staining for EpCAM (b) of cultured CTCs.




FISH (Fluorescent in situ hybridization) method for the detection of ALK rearrangement



Isolated red pattern predominant

Frequent separate red and green pattern with occasional isolated red

FISH (Fluorescent in situ hybridization) method for the detection of ALK rearrangement in cultured CTCs from lung cancer patient



Separate red and green pattern

## Conclusion

In our study, we showed that all FISH results in primary tumors were corresponding with real-time PCR results in cultured CTCs. Use of the cultured CTCs for molecular analysis has a merit of non-invasiveness and it could be easily repeated at different time-points during treatment to guide therapeutic decisions in a patient's treatment course. For the successful application of this strategy to clinical practice, CTC culture conditions will be further optimized. In addition, further confirmative characterization methods, such as different cell marker staining and molecular profiling, should be developed for precise identification of cultured CTCs.

## References

- [1] Cristofanilli M *et al.* Circulating tumor cells, disease progression, and survival in metastatic breast cancer. *N Engl J Med* (2004) 351:781-91.
- [2] Cohen SJ *et al.* Relationship of circulating tumor cells to tumor response, progression-free survival, and overall survival in patients with metastatic colorectal cancer. *J Clin Oncol* (2008) 26:3213-21.
- [3] Cristofanilli M *et al.* Circulating tumor cells: A novel prognostic factor for newly diagnosed metastatic breast cancer. *J Clin Oncol* (2005) 23:1420-30.
- [4] de Bono JS *et al.* Circulating tumor cells predict survival benefit from treatment in metastatic castration-resistant prostate cancer. *Clin Cancer Res* (2008) 14:6302-9.
- [5] Yu M *et al.* Circulating tumor cells: approaches to isolation and characterization. *J Cell Biol* (2011) 192:373-82.
- [6] Zer A and Leigh N. Promising targets and current clinical trials in metastatic non-squamous NSCLC. *Front Oncol* (2014) doi: 10.3389/fonc.2014.00329.
- [7] Faugeroux V *et al.* Clinical utility of circulating tumor cells in ALK-positive non-small-cell lung cancer. *Front Oncol* (2014) doi: 10.3389/fonc.2014.00281.